

Experiment 9: Water Purification: Hardness Estimation by EDTA method and its Removal using Ion-exchange Resin


Experimental values for the following students:

Titration-I: Standardization of EDTA

Pipette out 20 mL of the standard hard water containing 1mg/mL of CaCO_3 (1000 ppm) into a clean conical flask. Add one test tube full of ammonia buffer ($\text{NH}_4\text{OH} - \text{NH}_4\text{Cl}$) solution to maintain the pH around 10. Add three drops of Eriochrome Black – T (EBT) indicator and titrate it against the given EDTA solution taken in the burette. The end point is change of colour from wine red to steel blue. Repeat the titration for concordant titer values. Let ' V_1 ' be the volume of EDTA consumed.

S. No.	Volume of standard hard water (mL)	Burette reading (mL)		Volume of EDTA (V_1 , mL)
		Initial	Final	
1	20	0	18.5	18.5
2	20	0	19.0	19.0
3	20	0	18.5	18.5
Concordant titer value				18.5

Calculation:



Handwritten calculation showing the standardization of EDTA:

$$\begin{aligned} &20 \text{ mL of given hard water consumes } 18.5 \text{ mL of EDTA} \\ &20 \text{ mg of } \text{CaCO}_3 \text{ requires } 18.5 \text{ mL of EDTA} \\ &\therefore 1 \text{ mL of EDTA requires} = \frac{20}{18.5} \text{ mg } \text{CaCO}_3 \\ &= 1.081 \end{aligned}$$

This relation will be used in other two titrations

Titration-II: Estimation of total hardness of hard water sample

Pipette out 20 mL of the given sample of hard water into a clean conical flask. Add one test tube full of ammonia buffer ($\text{NH}_4\text{OH} - \text{NH}_4\text{Cl}$) solution and three drops of Eriochrome Black-T (EBT) indicator. Titrate this mixture against standardized EDTA solution taken in the burette. The end point is the change of color from wine red to steel blue. Repeat the titration for concordant titer value. Let ' V_2 ' be the volume of EDTA consumed.

S. No.	Volume of sample hard water (mL)	Burette reading (mL)		Volume of EDTA (V ₂ , mL)
		Initial	Final	
1	20	0	15.0	15.0
2	20	0	15.0	15.0
3	20	0	14.5	14.5
Concordant titer value				15.0

Calculation:

From titration -1 we have the following relation:

1 ml of EDTA requires = $20/v_1$ mg CaCO₃ for complexation

From titration -2

20 ml of sample hard water contains = v_2 ml of EDTA

= $v_2 \times 20/v_1$ mg of CaCO₃

\therefore 1000 ml of hard water sample consumes = $v_2 \times \frac{20}{v_1} \times \frac{1000}{20}$

$$= \frac{v_2}{v_1} \times 1000 \text{ ppm}$$

\therefore Total hardness of water sample = $\frac{15}{18.5} \times 1000 \text{ ppm}$

$$= 810.8 \text{ ppm (x)}$$

Titration-3: Removal of hardness using ion exchange method

Arrange the ion exchange column on to a burette stand and place a clean funnel on top of the column. Pour the hard water sample (around 40 to 50 mL) remaining after the completion of Titration – 2 through the funnel and into the ion exchange column. Place a clean beaker under the column and collect the water passing through the column over a period of 10 minutes. Adjust the valve of the column to match the duration of outflow.

From the water collected through the column, pipette out 20 mL into a clean conical flask and repeat the EDTA titration as carried out above. Note down the volume of EDTA consumed as 'V₃'.

Calculation:

From titration 1, we have the following relation

1 mL of EDTA requires = 1.081 mg of CaCO_3 for complexation

∴ 20 mL after softening through the softening column = 2.6 mL of EDTA

$$= 2.6 \times \frac{20}{18.5}$$

$$= 2.811 \text{ mg of } \text{CaCO}_3$$

$$\therefore 1000 \text{ mL of softened water} = \frac{V_3}{V_1} \times \frac{20}{20} \times \frac{1000}{20}$$

$$= \frac{V_3}{V_1} \times 1000 \text{ ppm}$$

$$= \frac{2.6}{18.5} \times 1000$$

$$= 140.5 \text{ ppm (4)}$$

S. No.	Volume of sample hard water (mL)	Burette reading (mL)		Volume of EDTA (V ₃ , mL)
		Initial	Final	
1	20	0	3	3
2	20	0	2.6	2.6
3	20	0	2.6	2.6
Concordant titer value				2.6

Result:

Total hardness of water sample = 810.8 ppm

Residual hardness in the water sample = 140.5 ppm

Hardness removed through column = 670.3 ppm