Abhinav Dinesh Srivatsa

Chemistry Lab

### **FAT Exam**

### Question Paper - Data Set 10

#### Part I

## Part-1

- 1. Principle of this experiment is that the rate of a reaction is dependant on a rate factor and proportional to the concentration of the reactants
- 2. The order of a reaction is the relationship the rate of a reaction and the concentration of it's reactants
- 3. The reagents used are: ethyl aretate, McI, water, Na On, ice when and the indicator phenolphthalien.
- 4. The ice when are used to immediately woll the solution and slow down the rate of reaction to take better measurements.
- 5. The application of this experiment is to find the titre values and see now temperature affects the rate of reactions relation to the concentration of ethylacetate.

### Part II

# Part 11

		Vol. of			
ζ. No.	Time (win)	Naon (unl)	(mr) ND -N+	109(4,)	$K = \frac{2.305}{+} \log \left( \frac{v_{\infty} - v_{\infty}}{v_{1}} \right)$
,	0	27.5	24.7	1. 393	0
2	10	27.8	24.2	1.384	0.00205
3	20	28.3	23.7	1.375	0.00207
4	30	28.5	23.5	1.371	0.00167
5	40	29.1	22.9	1.360	0.00 189
6	50	29.5	22.5	1.352	0.00187
7	10	520	-	-	-

Average K = 0.00205 + 0.00207 + 0.00167 + 0.00189 +0.00187

5

= 0.00191 min-1

Result

Cal whate value = 0.00191 min-1

Graphical value = 0.00189 min-1

Molewlarity = 2 (bimolewlar)

order of reaction = first order (psue do 1st order)

