


Fig 4.1. wave paths diagram.



$D$  = distance from the grating to the screen

$d$  = spacing between every two lines

If there are  $n$  lines per mm of grating, then

$$d = \frac{1}{n}$$

The diffraction grating formula for the principal maxima is:  $d \sin \theta = n \lambda$ , where  $n$  is the order of diffraction ( $= 1, 2, 3 \dots$ ) and  $\theta$ , the angle of diffraction

$$\lambda = \frac{d \sin \theta}{n}$$

Monochromator in  
Sophisticated Instruments

Apparatus Available:

Laser source, grating, scale with measurements

Student Learning Objectives:

To determine the wavelength of the given laser source using transmission diffraction grating method.

Observations:

Number of lines per meter on the grating is 98425.2.

Readings recorded in table 4.1.

Calculations:

$$\begin{aligned}\text{Mean wavelength} &= \frac{681 + 665 + 674}{3} \\ &= 673 \text{ nm}\end{aligned}$$

Result:

The wavelength of the laser source is found to be 673 nm.




Fig 4.2. Experimental Setup

$n$	$S$ (cm)	$2L$	$L$ (cm)	$\tan \theta = L/S$	$\theta = \tan^{-1}(L/S)$	$\sin \theta$	Mean	$\lambda$ (nm)
1	30	4	2	0.067	3.83	0.067	0.067	681
	40	5.4	2.7	0.068	3.89	0.068		
	50	6.6	3.3	0.066	3.78	0.067		
2	30	7.8	3.9	0.13	7.41	0.129	0.131	665
	40	10.6	5.3	0.133	7.58	0.132		
	50	13.3	6.7	0.134	7.63	0.133		
3	30	11.7	5.9	0.197	11.41	0.199	0.199	674
	40	16.2	8.1	0.203	11.48	0.199		
	50	20.2	10.1	0.202	11.42	0.198		

Table 4.1. Observations and calculations