

Upload the answers.
OBSERVATION AND CALCULATIONS


Potentiometric Titration-I:

Burette: KMnO_4 solution (0.05 N)

Beaker: 20 mL of steel solution containing Fe(II) + 20 mL (one test tube) of dil. H_2SO_4

Electrodes: Indicator electrode (Pt) to red terminal and SCE to black terminal

S. No.	Volume of KMnO_4 (mL)	EMF (volts)
1	0	0.05
2	1	0.06
3	2	0.07
4	3	0.08
5	4	0.16
6	5	0.16
7	6	0.17
8	7	0.18
9	8	0.19
10	9	0.19



Iron in carbon steel by potentiometry

Fig. 3: Plot of EMF vs Volume of KMnO_4 added (mL)

Potentiometric Titration-II:

Burette: KMnO_4 solution (0.05 N)

Beaker: 20 mL of steel solution containing Fe(II) + 20 mL (one test tube) of dil. H_2SO_4

Electrodes: Indicator electrode (Pt) to red terminal and SCE to black terminal

Sl. No.	Vol. of KMnO_4 (mL)	EMF (Volt)	ΔE (Volt)	ΔV (mL)	$\Delta E/\Delta V$ (Volt/mL)	Average Volume (mL)
1	3.0	0.08				
2	3.2	0.09	0.01	0.2 ml	0.05	3.1 ml
3	3.4	0.095	0.005	0.2 ml	0.025	3.3 ml
4	3.6	0.15	0.055	0.2	0.275	3.5
5	3.8	0.16	0.01	0.2	0.05	3.7
6	4.0	0.161	0.001	0.2	0.005	3.9
7	4.2	0.163	0.002	0.2	0.01	4.1
8	4.4	0.165	0.002	0.2	0.01	4.3
9						
10						




Fig. 4: Plot of $\Delta E/\Delta V$ vs Average volume of KMnO_4 added.

Calculation:

Iron in carbon steel by potentiometry

$$(N \times V) \text{ of steel sample solution} = (N \times V) \text{ of KMnO}_4$$

$$\begin{aligned} N \text{ of steel sample solution} &= \frac{0.05 \text{ N} \times \text{Volume of KMnO}_4 \text{ from Plot-2}}{20 \text{ mL of steel sample}} \\ &= \underline{0.00875} \text{ N} \end{aligned}$$

Amount of Fe present in 1 L of sample solution = Normality of steel sample x At. wt. of Fe (55.85)

$$\begin{aligned} \text{Amount of Fe present in given (100 ml) sample solution} &= \frac{\text{Normality of steel sample} \times 55.85 \times 100}{1000} \\ &= \underline{0.049} \text{ grams in 100 mL} \end{aligned}$$

Result: The amount of Iron present in given steel sample is found to be = 0.0098 grams.

$\Delta E / \Delta V$ (volt/mL)

Scale

X-axis \rightarrow $\text{vol} = 0.1 \text{ mL}$

Y-axis \rightarrow $\text{Zem} = 0.04 \text{ volt/mL}$

0.32

0.28

0.24

0.20

0.16

0.12

0.08

0.04

0

3.1

3.2

3.3

3.4

3.5

3.6

3.7

3.8

3.9

4.0

4.1

4.2

4.3

4.4

4.5

Average volume of KMnO_4 (mL)

